



# C-QUR™ Mesh

Available in the UK from  
Pierson Surgical Ltd  
sales@piersonsurgical.com  
01225 766632  
www.piersonsurgical.com



Omega 3 fatty acid (O3FA) coated polypropylene mesh designed for open and laparoscopic hernia repair when a tissue separating layer is needed.

C-QUR™ Mesh combines Atrium's ProLite™ polypropylene mesh with a tissue separating layer of Omega 3 fatty acid (O3FA). C-QUR Mesh is designed for intra-abdominal placement.

- Strong knit construction<sup>1</sup>
- See-thru clarity
- O3FA coating has been shown to minimize tissue attachment to the mesh<sup>2</sup>
- Ability to trim allows surgeon the option to customize the mesh

#### Competitive Device

- Davol™/Bard®: Composix™ E/X & LP
- Davol™/Bard®: Sepramesh™ IP
- Davol™/Bard®: Ventralight™ ST
- Ethicon: Proceed™
- Ethicon: PhysioMesh™
- Covidien: Parietex™ Composite & Parietex™ Composite Optimized
- Gore™: DualMesh




1. Deeken CR, Abdo MS, Frisella MM, Matthews BD. Physicomechanical evaluation of polypropylene, polyester, and polytetrafluoroethylene meshes for inguinal hernia repair. *J Am Coll Surg*. 2011 Jan;212(1):68-79.

2. Pierce, R.A. et al. "120-Day Comparative Analysis of Adhesion Grade and Quantity, Mesh Contraction, and Tissue Response to a Novel Omega-3 Fatty Acid Bioabsorbable Barrier Macroporous Mesh After Intraperitoneal Placement." *Surg Innov* (2009) Volume 16, Number 1, 46-54

## PRODUCT SIZES AND CODES


### Round

Order Qty.	Part No.	Size (cm)	Size (in)	Box Qty.
	31505	9	3.54	1 each
	31506	13	5.12	1 each
	31507	15	5.90	1 each




### Square

Order Qty.	Part No.	Size (cm)	Size (in)	Box Qty.
	31543	8.9 x 8.9	3.5 x 3.5	1 each
	31544	12.5 x 12.5	5 x 5	1 each




### Oblong/Oval

Order Qty.	Part No.	Size (cm)	Size (in)	Box Qty.
	31526	7.5 x 15	3 x 6	1 each
	31528	10 x 15	4 x 6	1 each
	31527	11 x 14	4.3 x 5.5	1 each
	31531	13.5 x 17.8	5.4 x 7	1 each
	31533	15 x 20.3	6 x 8	1 each
	31534	15 x 23	6 x 9	1 each




### Mid-Line/Oval

Order Qty.	Part No.	Size (cm)	Size (in)	Box Qty.
	31538	16 x 31.2	6.3 x 12.3	1 each



### X-Large Oblong/Oval

Order Qty.	Part No.	Size (cm)	Size (in)	Box Qty.
	31535	17.8 x 23	7 x 9	1 each
	31536	20.3 x 25.4	8 x 10	1 each
	31537	20.3 x 30.5	8 x 12	1 each
	31540	25.4 x 35.5	10 x 14	1 each
	31542	30.5 x 45.7	12 x 18	1 each



Available in the UK from  
 Pierson Surgical Ltd  
 sales@piersonsurgical.com  
 01225 766632  
 www.piersonsurgical.com



© Atrium Medical Corporation 2013. All rights reserved. 09/13 Part # 0670.

Atrium Medical Corporation • 5 Wentworth Drive, Hudson, NH 03051 • Phone: 603-880-1433 Fax 603-880-6718.

CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. Refer to Instructions for Use for current indications, warnings, contraindications, and precautions. Atrium, C-QUR and ProLite are trademarks of Atrium Medical Corporation, a MAQUET GETINGE GROUP company.

All third party trademarks used herein are property of their respective owners.

Atrium products are protected by international and U.S. Patent(s). For patent information visit: <http://patents.maquet.com>