

FEATURES



- No memory
- Increased strength
- Specific extrusion process
- Optimal smooth

BENEFITS

- ▶ Easy to handle
- ▶ Safety
- ▶ Perfect slide
- ▶ Reduced tissue trauma

USES

- CABG procedure
- Vascular surgery



PACKAGING

- Coropak®, 4-Pak® for minimal memory
- Specific Teflon® pledgets




RANGE

- USP 9/0 to 0
- Needle length : 6,4 mm to 36 mm
- Thread length : 45 cm to 120 cm

STUDIES

- Thinner and Stronger than competitors

- Special high quality extrusion process



- Corolene®

USP 8/0



- Competitor

Source : Service R&D / R&D department

Strand breaking strength (N)

Comparison Chart - Péters Surgical Non Absorbable Sutures

	PREMIO®	CO R O L E N E®	CARDIONYLL®	CARDIOFLON®	CARDIOXYL®
Structure	Monofilament	Monofilament	Monofilament	Braided	Braided
Material	PVDF	Polypropylene	Polyamide	Polyester impregnated with Teflon®	Silicone Coated Polyester
Key Features	<ul style="list-style-type: none"> • Excellent long term stability • Minimal memory • Excellent knot strength • Optional Pledgets 	<ul style="list-style-type: none"> • Low friction, superior glide , particularly suited to the “parachute” technique • Very low memory easy to handle • Minimal tissue trauma • High strength 	<ul style="list-style-type: none"> • Minimal memory • Very supple • Guaranteed knot strength • Optimal elasticity - reduces the number of knots required • Optional Pledgets 	<ul style="list-style-type: none"> • Very high strength • The Teflon induction gives very high knot security • Superior suppleness • Optional Pledgets 	<ul style="list-style-type: none"> • Silicone coating • Slides through tissue with ease • Simplifies knot adjustability and positioning
Colour	Blue	Blue	Blue	Green or Green & White	Green & White
USP Range	9/0 – 2/0	8/0 – 2/0	6/0 – 1	5/0 - 5	5/0 – 2/0
Surgical use	Neonatal and Paediatric Surgery Vascular Surgery Cardiac Surgery	Cardiac and Vascular Surgery Coronary Surgery	Cardiac Surgery Mitral and Aortic Valve repair Cardiac cavity closures	Cardiac Valve Repair Cardiac valve and ring placement Securing CPB Cannula Thoracic Drain placement	Cardiac valve and ring placement Securing CPB Cannula Patch implantation

Pierson
Surgical Ltd.



Distributed in the UK by Pierson Surgical

Email: sales@piersonsurgical.com

Phone: 01225 766632 Fax: 07092 315510

www.piersonsurgical.com

Pierson Surgical Ltd, North Bradley House, North Bradley, Trowbridge, BA14 0TA